

User Manual for Retinal Disease Detection System

Comprehensive Guide for Patients and Doctors

Table of Contents

1. Introduction
 - 1.1. System Overview
 - 1.2. Key Features
 - 1.3. Who Should Use This System?
2. Getting Started
 - 2.1. System Requirements
 - 2.2. Creating an Account (Registration)
 - 2.3. Logging In
 - 2.3.1. Password Recovery
3. Uploading and Analyzing Fundus Images
 - 3.1. Image Requirements & Best Practices
 - 3.2. Step-by-Step Upload Process
4. Recommendations
5. Viewing History & Reports
6. Feedback & Support
 - 6.1. Overall Feedback
 - 6.2. Individual Feedback
 - 6.3. Contacting Admin
7. Troubleshooting & FAQs
8. Glossary of Terms
9. Privacy & Security

1. Introduction

System Overview

This web-based medical imaging system helps detect diabetic retinopathy and glaucoma by analyzing fundus images of the retina. It provides instant results, recommendations, and a secure history log for tracking eye health.

Key Features

- ✓ User-friendly interface for patients & doctors
- ✓ Automated disease detection with AI analysis
- ✓ Recommendations based on results
- ✓ Secure image storage & history tracking
- ✓ Direct feedback & admin support

Who Should Use This System?

- Patients monitoring their eye health
- Doctors seeking a quick diagnostic aid
- Clinics integrating AI-based screening

2. Getting Started

System Requirements

- A modern browser (Chrome, Firefox, Edge)
- Stable internet connection
- High-quality fundus image (JPEG/PNG, min. 1024x1024 pixels)

Creating an Account (Registration)

1. Click "Sign Up" on the homepage/ you can also sign up directly with google.
2. Fill in:
 - Name
 - Email
 - Password (8+ characters, including numbers & symbols)
3. Verify your email via the confirmation link sent.

Tip: Use a strong password and keep it secure.

Logging In

1. Visit the login page.
2. Enter your email and password / Sign in via google.
3. Click "Login".

Password Recovery? Click "Forgot your Password" and follow the email instructions.

3. Uploading and Analyzing Fundus Images

Image Requirements & Best Practices

- ✓ High resolution (1024x1024px minimum)
- ✓ Clear focus (no blur or reflections)
- ✓ Well-lit retina (avoid shadows)
- ✗ Avoid low-quality, cropped, or tilted images.

Step-by-Step Upload Process

1. After logging in, click "Choose Image".
2. Drag & drop or browse your fundus image.
3. The system will adjust image quality.
4. If approved, click "Upload Images".
5. If rejected, retake/upload a clearer image.

4. Recommendations

Based on your results, the system will suggest:

- ✓ Next steps (e.g., see a retina specialist or glaucoma workup if needed)
- ✓ Lifestyle tips (control blood sugar, blood pressure, avoid smoking)
- ✓ Recommended follow-up interval for future screenings

5. Viewing History & Reports

Go to "History" in the dashboard.

You can view your all past scans.

6. Feedback & Support

Overall Feedback:

Click "Feedback" in the menu.

Rate your experience and describe any issues.

Submit – the admin will review your comments.

Individual Feedback: You can provide feedback on each of your results on the home page in the results section or on the history page.

Contacting Admin:

Please submit your query in writing, providing clear details of your issue or question, and our support team will respond promptly.

8. Glossary of Terms

Fundus Image: A photograph of the retina.

Diabetic Retinopathy: Diabetes-related eye damage.

Glaucoma: A condition damaging the optic nerve.

9. Privacy & Security

- Your data is encrypted.
- Images are stored securely.
- Only you can access reports.

- This manual ensures smooth navigation through the retinal disease detection system. For further assistance, contact our support team.

Stay proactive about your eye health!

Last Updated: 1/07/2025

© 2024 Retinal Disease Detection System. All rights reserved.